

Item No.	DMS-T2	
Capacity	2.0L	
Product Size	Length	19cm
	Width	11cm
	Height	16.5cm
Carton Size	40.5x27x45cm	
Quantity in one Carton	30pcs	
G.W. / N.W.	5.5KGS / 4.5KGS	

2.0L T-SERIES SHARPS CONTAINER

The 2.0L Sharps Container T-Series is a popular small-size Biohazard container for safe storage and easy manipulation (2.0L capacity). Case of 30.

Regulatory Approvals: ISO 9001:2008 Certified; Products are BS7320, UN and FDA approved.




PRODUCT DESCRIPTION

The DMS-T2 Sharps Container 2.0L is a great small T-series Biohazard container with a secure perimeter lock for safe transportation and counterbalanced door for hand-free disposal.

Features:

- Most popular quadrate shape
- Leak proof and puncture resistant
- Wide and clear opening for easy disposal
- Counterbalanced door for hand-free disposal
- Ergonomically-designed handle
- Temporary closure & definite closure
- Practical rotating/sliding temporary lid
- Secure perimeter lock for safe usage and transportation
- Easy manipulation
- Comprehensive label details
- Safe to human health and environment



Item No.	DMS-T2	
Capacity	2.0L	
Product Size	Length	19cm
	Width	11cm
	Height	16.5cm
Carton Size	40.5x27x45cm	
Quantity in one Carton	30pcs	
G.W. / N.W.	5.5KGS / 4.5KGS	

Dailymag is ISO 9001:2008 Certified, and our products are BS7320, UN and FDA approved. We have more than 100 patents, including more than 20 for sharps containers. We are China’s leading manufacturer and exporter of sharps containers and Sharps Boxes.<http://www.sharpcontainer.com/>10 Business Days. Recommended minimum order of \$2,000 in combined sharps container products to qualify for freight shipping.

- (x30) DMS-T2 Sharps Container 2.0L

What is Sharps Container?

Sharps container protects workers from the potential of being struck by loose and exposed sharps and helps to provide convenient and affordable disposal services to the public

What are sharps?


Sharps include needles, IV tubing with needles attached, scalpel blades, lancets, glass tubes and slides, and syringes. According to international standard, approved sharps container are rigid, leak-proof, puncture resistant, sealed and clearly marked with the bio-hazard symbol.

Why use the Sharps containers?



GLOBAL HEALTH

2.0L T-SERIES
SHARPS CONTAINER



The first priority in Sharps container's development is safety ... for the patients, doctors and for personal disposal. All sharps containers are designed for safe use, effectively blocking hand access from discarded sharps. Units snap together, permanently locking with tight seams between the base and covers. Walls are puncture resistant, meeting or exceeding international Standards.

- For optimum efficiency, a variety of sharps disposal units are offered to clients reducing overall waste.

Rigid polypropylene construction with handle

- Safe to incinerate and environmentally friendly
- Puncture and impact resistant
- Fluid leak resistant
- Translucent top helps prevent overfilling & injury
- Tough and durable with strong handles
- Secure final lock prevents re-opening
- Needle removable notches at top
- Inter-stackable minimizes storage & ease transportation
- Mountable also on rail, IV pole, bed and trolley