



AJ-06B MULTI-PURPOSE DELIVERY BED

This table can be used for pregnant examine, maternity childbearing, dystocia surgical and curettage gynecologic surgery.

PRODUCT DESCRIPTION

This table can be used for pregnant examine, maternity childbearing, dystocia surgical and curettage gynecologic surgery. This table is raised and lowered by the hydraulic pressure, the front and back lean, backboard fold is the worm, gear, rack structure, Various actions can be adjusted in the right scope and easy to operate, The leg board can be disassembled, All the outside cover and ordure basin are stainless steel , it is specious and easy to clean. **Specifications**

- Length: 1330□1880mm
- Width: 600mm
- Height: 680□930mm
- Adjust ≥ 240 mm
- Reversed Trendelenburg $\geq 10^\circ$
- Trendelenburg $\geq 20^\circ$
- Back Section Raised From the Horizontal $\geq 45^\circ$
- Back Section lowered From the Horizontal $\geq 10^\circ$