

AJ-2301 INFANT INCUBATOR

Simultaneously air temperature and set temperature display.

PRODUCT DESCRIPTION

Features

- Simultaneously air temperature and set temperature display.
- Alarm functions for safety includes
- Over temperature deviation
- Temperature probe failure
- $>37^{\circ}\text{C}$ temperature setting
- A keypad locks helps to prevent inadvertent changes to settings.
- Trough for pulling and pushing add the water and wash is convenient.
- Heater is electrically cut off if air /skin temperature inside incubator exceed 40°C .
- The infant bed inclination angle is adjustable.


Technical specifications

- Electrical requirement: $220\text{V}\pm 10\%/50\text{HZ}\pm 2\%$
- Power supply : $\leq 420\text{VA}$
- Air Temperature control range $\square 25^{\circ}\text{C}\sim 38^{\circ}\text{C}$
- Temperature fluctuation: $\leq 0.5^{\circ}\text{C}$
- Uniformity of mattress temperature: $\leq 0.8^{\circ}\text{C}$
- Warm-up time: $\leq 45\text{min}$
- Internal noise level: $\leq 55\text{dB}$



GLOBAL HEALTH

AJ-2301 INFANT
INCUBATOR



Operation condition

- Environment temperature: 25°C~38°C
- Environment relative humidity: 30%~75%RH
- Environment air velocity of flow: <0.3m/s

Package

- Gross weight: 110kg
- Net weight: 80kg
- Measure: 1100*780*860mm