



## AJ-3000CT 12.1-INCH TOUCH SCREEN PATIENT MONITOR

With a variety of optional extras, including arrhythmia and ST analysis, EtCO<sub>2</sub> measurement and double IBP, AJ-3000CT touch screen multi-parameter patient monitor is configurable to suit your requirements.

### PRODUCT DESCRIPTION

Purchase only the advanced features that you need and nothing more. With NIBP recall and extensive data storage capability, AJ-3000CT touch screen multi-parameter patient monitor is particularly well-suited for use in peri-anesthesia and critical care environments. Rich in functions, conveniently adaptable and pleasingly affordable, AJ-3000CT touch screen multi-parameter patient monitor may be just the monitor for your department. **Features & Advantages**

- 6 Standard parameters: ECG, RESP, NIBP, SPO<sub>2</sub>, 2-TEMP, PR/HR
- 12.1"high resolution color touch screen display
- Real time S-T segment analysis, pace-maker detection
- Indication of Incorrect manipulation and malfunction analysis
- Patient info input management function
- Multi-lead ECG waveforms display in phase
- Large volume of tabular and graphic trends information storage and easy to recall
- Capture dynamic waveforms
- Efficient resistance to interference of defibrillator and HF knife
- UP to 2 hours working capacity of built-in rechargeable battery
- Unique accessory management cabinet
- Networking capacity
- Option: ETCO<sub>2</sub>, IBP, and 2-IBP