



AJ-3002 MULTI-PURPOSE O.T. TABLE (SIDE CONTROLLED)



This table is suitable for surgical procedures of the head, neck, thorax and abdomen, perineum and extremities as well as for gynecology and orthopedics operation.

PRODUCT DESCRIPTION

This table is suitable for surgical procedures of the head, neck, thorax and abdomen, perineum and extremities as well as for gynecology and orthopedics operation. The table's elevation, lowering, reversed trendelenburg, trendelenburg, lateral tilt can be all automatically controlled. Leg section can be rotated, disassembled and adjusted easily. It is most suitable for the urological operation. The table top made of high strength glass, can move backwards and forwards. X-Ray examination can be carried out on it, It can be raised or lowered hydraulically by means of an oil pump. Various positions can be controlled from both sides.

Chassis inside an outside covers are stainless steel. **Specifications**

- Length: 2100mm
- Width: 480mm
- Max Height: 1040mm
- Height: 800mm
- Adjust ≥ 240 mm
- Reversed Trendelenburg $\geq 20^\circ$
- Left $\geq 18^\circ$
- Right $\geq 18^\circ$
- Trendelenburg $\geq 20^\circ$
- Head Section lowered from the Horizontal $\geq 90^\circ$
- Head Section Raised from the Horizontal $\geq 45^\circ$
- Backward and Forward Movement ≥ 400 mm
- Leg section Opened $\geq 90^\circ$
- Back Section Raised From the Horizontal $\geq 55^\circ$
- Leg section Lowered from the Horizontal $\geq 90^\circ$
- Kidney Bridge Elevation ≥ 120 mm