



AJ-50BY 50MA MOBILE X-RAY MACHINE

This X-ray apparatus features full wave rectification and modular X-ray tube head.

PRODUCT DESCRIPTION

This X-ray apparatus features full wave rectification and modular X-ray tube head. The frame of the unit, of cantilever construction, is compact in design easy of movement. The machine can be used for radiography in hospital wards or operating rooms. **Specifications**

- Power requirements
 - Voltage:180-240V Frequency:50Hz
 - Internal resistance: $\leq 1\Omega$ Rating:5k VA
 - Current: Maximum 15A for radiography
- Time: 0.063 to 5.0 steps electronic
- Maximum height from X-ray tube focus to floor >1720mm
- Minimum height from X-ray tube focus to floor <500mm
- Collimator: Maximum film size at 1000mm
- Focal distance: 400×400mm
- Maximum remote control distance: 6m
- X-ray tube focus: 2.6
- Unit moving force: $\leq 250N$
- Net weight of the unit:95kg