



AJ-H SLIDE DRYER

New heating material is heat conduct fast, durable, energy-saving.

PRODUCT DESCRIPTION

Features

- New heating material is Heat conduct fast, durable, energy-saving
- Display presetting temperature & current temperature
- Self memory of set temperature after running
- Corrosion-resistance of the bath

Technical Specifications

- Range of setting temperature: ambient to 90°C (constant adjusting)
- Temperature tolerance: + 1°C
- Allowance ambient temperature: 0 - 40°C
- Working voltage: AC220V/50Hz, AC110V/60Hz
- Power: 200 VA
- Dimension: 340 x 340 x 120 mm
- Weight: 6 kg