



GLOBAL HEALTH

AJ-S900D HIGH
FREQUENCY
ELECTROSURGICAL
MACHINE



AJ-S900D HIGH FREQUENCY ELECTROSURGICAL MACHINE

Range of Use: May use to various surgical operations of cutting or/and coagulation to organism tissue.

PRODUCT DESCRIPTION



Specifications Range of Use: May use to various surgical operations of cutting or/and coagulation to organism tissue. Working Modes: Pure cut, Blend 1,Blend 2, Blend 3, Desiccate, Fulgurate, Spray.

Working Modes	Pure cut or cut	5□300W(may provide type 5□400W)
	Blend 1	5□240W
	Blend 2	5□200W
	Blend 3	5□150W
	Desiccate	5□200W
	Fulgurate or Coag.	5□150W
	Spray	5□100W
	Bipolar	-
Accessories	Active Electrode	5
	Footswitch	1
	Neutral Electrode	1
	Bipolar Electrode	-
Transport Data	Disposable Electrode	Pick out and buy
	Package Dimensions	45×26.5×47cm
	Gross Weight	11kg
	Voltage	220V ±22V,50Hz/60Hz
Technical Data	Working Frequency	512kHz ±5%
	Variety of Equipment	class 1,type CF.Ordinary Equipment
	Working Fashion	Loading intermittently, Functioning continuously.
	Leskage Currents To Ground	Less than 0.5mA
	Leskage Currents Through Patient	Less than 0.01mA