




GLOBAL HEALTH

AJ-TG12M
MICROHEMATOCRIT
CENTRIFUGE



AJ-TG12M MICROHEMATOCRIT CENTRIFUGE

This product can be used to determine the specific product value of blood cell and separate the micro-blood from micro-solution.

PRODUCT DESCRIPTION

Features This product can be used to determine the specific product value of blood cell and separate the micro-blood from micro-solution. **Technical Parameters**

- Max speed: 12000r/min
- Maximum Relative Displacement: 14800xg
- Rotor capacity: 24 pieces of capillaries
- Time range 0-99min
- Motor: DC brushless motor, Microprocessor control
- Power supply: AC220V 50Hz 2A
- Noise: ≤62dB(A)
- Dimension: 240×260×250mm(L×W×H)
- Weight: 12 kg