



AJ-TS3A AUTOMATIC TISSUE PROCESSOR

The intelligent design gives the ability to deal with abnormal situation automatically.

PRODUCT DESCRIPTION

Features

- The intelligent design gives the ability to deal with abnormal situation automatically
- The imported blue LCD is clear & stable with several control system
- Recirculated atmosphere depurative system that absorbing unhealthy volatilization gas
- Anhydrous heating element with high accuracy
- UPS build-up, power backup for continuous working after power off
-

With key parts imported, it runs steadily without noise and comply with theory of ergonomics

Technical Specifications

- Numbers of tanks: 12 (9 liquid tanks , 3 paraffin tanks)
- Process time: 0~99 hrs for first tank , 0~24 hrs for other tanks
- Paraffin tank temperature: ambient temperature to 99°C, double protection
- Heating Method: Waterless inner heating
- Capacity : 1000ml
- Temperature tolerance : $\pm 1^{\circ}\text{C}$
- Agitate times: 2/min
- Power off protection (it can working ~6 hrs after power off)
- Room temperature: 0~40°C
- Working voltage: AC220V/50Hz, AC110V/60Hz
- Power Draw: 800 W
- Dimension: 1000 x 420 x 450 mm