



AJ004 ICU BED / ELECTRIC BED MULTI FUNCTION

Integrated stamping bed-tops and frame are steel made with epoxy polyester powder coated & baked finish, resistant to chipping and scratching.

PRODUCT DESCRIPTION

Specifications

- Electric bed with I.V. Pole, powder coated steel frame
- The removable headboards are made of ABS plastic
- Integrated stamping Bed-Tops and frame are steel made with epoxy polyester powder coated & baked finish, resistant to chipping and scratching
- 6-stuff aluminum alloy guardrail protect patient even back and keen rest have been raised
- Mounted on four 125mm diameter with brake
- The back panel can be lifted up from 0° to 80° ($\pm 5^\circ$), the foot panel can be lifted up from 0° to 45° ($\pm 5^\circ$), height adjustment 250mm, trendelenburg and reversed trendelenburg 12°, 2120*980*500-750mm
- Matched with I.V. Pole
- 2, 3, 5 functions for your choice
- 2 Function: Back-rest lifting, Leg-rest lifting
- 3 Function: Back-rest lifting, Leg-rest lifting, Height adjustment
- 5 Function: Back-rest lifting, Leg-rest lifting, Height adjustment, Trendelenburg and Reversed Trendelenburg.