




GLOBAL HEALTH

AJ006 ICU BED /
ELECTRIC BED
MULTI FUNCTION



AJ006 ICU BED / ELECTRIC BED MULTI FUNCTION

Integrated stamping Bed frame is steel made with epoxy polyester powder coated and baked finish, resistant to chipping and scratching; Bed surface is ABS material.

PRODUCT DESCRIPTION

Specifications

- Electric bed with I.V. Pole and urine hook
- The removable head and foot boards are made of ABS
-

Integrated stamping Bed frame is steel made with epoxy polyester powder coated and baked finish, resistant to chipping and scratching; Bed surface is ABS material

- ABS 4-piece tuck-away guardrails protect patient even back and knee rest have been raised
- Mounted on 4*150mm castors, with central lockable system
- Matched with I.V. Pole
- 2, 3, 5 functions for your choice
- 2 Function: Back-rest lifting, Leg-rest lifting;
- 3 Function: Back-rest lifting, Leg-rest lifting, Height adjustment;
- 5 Function: Back-rest lifting, Leg-rest lifting, Height adjustment, Trendelenburg and Reversed Trendelenburg.

Functions:

- The back rest can be lifted up from 0° to 80° ($\pm 5^\circ$);
- The leg rest can be lifted up from 0° to 45° ($\pm 5^\circ$);
- Height adjustment 250mm;
- Trendelenburg and reversed trendelenburg 12°.
- Dimension: 2080 x 1020 x 500-750mm