



## DHP-050 50L INCUBATOR AIR OVEN

Specular stainless-steel chamber, semicircular arcs at corners for easy cleaning and the space between the shelves in the chamber is adjustable.

### PRODUCT DESCRIPTION

#### Features

- Specular stainless-steel chamber, semicircular arcs at corners for easy cleaning and the space between the shelves in the chamber is adjustable
- Microprocessor temperature controller ensures a precise and reliable control
- Forced warm-air circulation design ensures a uniform temperature in chamber
- Closed-silicon gasket and double glass door design make observation easily and clearly

#### General Specification

- Electrical requirements : 220V 50Hz
- Temperature range : as incubator: 5~ 80° C , as oven : 80~ 200° C
- Temperature control : 0.1° C
- Temperature stability : a.  $\pm 0.5^{\circ}$  C; b.  $\pm 0.5\%$
- Ambient temperature : 5~ 40° C
- Timing range : 1~ 9999min
- Input power : 850W
- W\* D\* H interior dimension ( mm ) : 420\* 350\* 350
- W\* D\* H overall dimension ( mm ) : 720\* 590\* 520
- Shelves per chamber : 2 pcs