



## DR181 SERIES DIGITAL HOLTER RECORDER, 12- LEAD



The DR181 has a 12-lead configuration in addition to its 3- and 5-lead versions. It is a high-quality research and clinical evaluation tool which records 12-lead / 9 channels of data without any use of arithmetically derived or reconstructed data from non-standard Holter lead hook-ups.

**Regulatory Approvals:** ISO, CE, FDA

### GALLERY IMAGES



### PRODUCT DESCRIPTION

The DR181 is the latest version of the reliable and popular DR180 series. The DR181 in a 12-lead configuration records 12-lead / 9 channels of data without shortcuts. The 12-lead system does not use arithmetically derived or reconstructed data from non-standard Holter lead hook-ups. It provides a high quality research and clinical evaluation tool. The DR181 12-lead settings can be adjusted to fit the 12-lead recording requirements, intervals and strip lengths are also adjustable.



DR181 SERIES  
DIGITAL HOLTER  
RECORDER, 12-  
LEAD



Features:

- Meets EC-11 and EC-38 standards;
- Ergonomic design and proven reliability;
- Removable digital storage memory (compact flash);
- Reimbursement for all applicable CPT (Current Procedural Terminology) codes for Holter;
- Shielded cables and Pacemaker detection standard;
- Foreign language interfaces;
- Industry-leading 3-year warranty;
- Cable options - 5-lead, 7-lead or Oxymetry (oxy plus 5-lead) cables;
- Maximize good quality recording cycles;
- Reliable digital recording;
- Pediatric (high sensitivity mode);
- Compatible with NorthEast Holter LX® Analysis software.

[Holter DR181 12-Lead Product Brochure](#)

[DR181 Holter Recorder Hook-up Guide](#)

[DR181 Operator's Manual \(Rev-F\)](#)

[Holter Warranty Terms](#)

- **12-Lead Patient Cable**
- **Recorder Pouch and Straps**
- **DR181 Battery Door Cover**
- **DR181 12-Lead Hook-up Kit**
- **DR181 Memory - SD Flashcard**
- **Nonin Sensor Finger Probe**
- **Nonin Flex-wrap pkg - 25**
- **Nonin Flex-Wrap (infant) pkg - 25**
- **Nonin Finger Clip**
- **Nonin Ear Clip**
- **Nonin Forehead Sensor**



NorthEast Monitoring is a medical electronics company that develops high quality, rugged products to address diagnostic, clinical and research needs in cardiology, respiratory, and pulmonary applications. <http://www.nemon.com/>

- (1) DR181 Holter, 12 channel digital recorder
- (2) 10-lead/12Channel patient cable
- (3) 10 wires for 10-wire/12 channel patient cable (wires only)
- (4) Pouch with shoulder and waist strap
- (5) 2GB SD Flash Memory Card
- (6) Users manual
- (7) Hook-up guide
- (8) Hook-up kit 7 business days

NorthEast Monitoring, Inc. warrants that the company's products meet the labeled specifications of the products and will be free from defects in materials and workmanship under normal consumer usage for the period(s) outlined below.

Warranty Period for:

- a) Holter Recorders, Event Recorders and combination Holter/Event Recorders is Three (3) years from the date of shipment from NorthEast Monitoring, Inc.
- b) DR180+ and DR181 Holter Recorder Shielded Patient Leads, including the DR180+ and DR181 Oximetry Shielded Patient Leads is Six (6) months from the date of shipment from NorthEast Monitoring, Inc.

Note: Due to the limited life of the oximetry sensor, this part of the cable system is not under warranty.

The warranty is void in cases of:

- a) Acquisition of Products through Unauthorized Channels (applies to products delivered outside the US only).
- b) Periodic maintenance, repair and replacement of parts due to normal wear and tear
- c) Defects or damage that result from improper operation, storage, misuse or abuse, accident or neglect, such as physical damage (cracks, scratches, etc.) to the surface of the product resulting from misuse; contact with liquid, water, rain, extreme humidity or heavy perspiration, sand, dirt or the like, extreme heat, or food; subjecting the Product or Accessory to abnormal usage or conditions.
- d) Use of Non-NorthEast Monitoring, Inc. Products and Accessories;
- e) Unauthorized Service or Modification;
- f) Products or Accessories with (a) serial numbers or date tags that have been removed, altered or obliterated;
- g) Software Embodied in Physical Media;
- h) Products that are Lost or Stolen.

If a product covered by this warranty is determined to be qualified under this warranty,

and the warranty claim is made within the warranty period, NorthEast Monitoring, Inc. will, at its discretion, repair or replace the defective part(s) free of charge.