




GLOBAL HEALTH

FLY-1102C 170L  
SHAKING  
INCUBATOR



## FLY-1102C 170L SHAKING INCUBATOR

### PRODUCT DESCRIPTION

- Rotary speed: 40~300 rpm
- Rotary speed precision:  $\pm 1$  rpm
- Vibration range:  $\phi 26$ mm
- Maximum capacity: 100mL \*48 or 250 mL\*25
- Standard configuration: 251 mL \*12
- Size of rocking plate: 500\*350 (mm)
- Timing scope: 0.999 hours
- Temp. control scope: 5~60°C
- Temp. control precision:  $\pm 0.1^{\circ}\text{C}$  (Homoeothermic)
- Temp. Evenness:  $\pm 1^{\circ}\text{C}$
- Display: LCD
- Number of rocking plates: 3 (accessories available)
- Shape and dimension: 702\*562\*1262 (mm)
- Packing size: 852\*742\*1432 (mm)
- Weight: 150kg
- Capacity: 602\*442\*652 (mm) (170L)
- Power: 436W
- Power supply: AC220V  $\pm 10\%$  50~60 Hz