





GLOBAL HEALTH

PQX-450H 450L
ARTIFICIAL CLIMATE
INCUBATOR



PQX-450H 450L ARTIFICIAL CLIMATE INCUBATOR

PRODUCT DESCRIPTION

General Specifications

- Temperature range : With Lighting 10~ 50° C, Without Lighting 4~ 50° C
- Temperature control : 0.1° C
- Humidity stability : $\pm 0.1^{\circ}$ C
- Electrical requirements : 220V 50Hz
- Ambient temperature : 5~ 35° C
- Continuous working time : No less than 180h
- Cubage : 450L
- Temperature uniformity : $\pm 2^{\circ}$ C
- Humidity range : 50~ 90% RH
- Humidity stability : $\pm 5\% \sim 7\%$ RH
- Lighting intensity : 0~ 25000LX
- Input power : 2100W
- W* D* H interior dimension (mm) : 700* 550* 1140
- W* D* H overall dimension (mm) : 950* 850* 1850
- Shelves per chamber : 4 pcs