



**Innovative**  
The LifeWrap is a low-cost, low-technology, first-aid device to treat hypovolemic shock



**Integrated**  
Can be used with other interventions, kept on during transportation procedures and surgery




**Practical**  
The suit is affordable, safe, and accessible-ideal for low resource settings


**Product Fact Sheet for the LifeWrap NASG**

Date:

<b>Description:</b>	<b>LifeWrap NASG</b> Non-Pneumatic Anti-shock Garment		
<b>Size</b>	<b>LARGE / GREEN</b>		
<b>Product Code:</b>	LWNASGLGE15		
<b>Manufacturing MOQ</b>	1,000 Units		
<b>FOB USD\$ FOB HK</b>			
<b>Below Production MOQ</b>	50 Units		
<b>FOB USD\$ HK</b>	POA		
<b>Bar Code #:</b>	0680569839954		
<b>Individual Packaging:</b>	Polybag		
<b>Ind Packed Measurements</b>	<b>L</b>	<b>W</b>	
<b>Inches:</b>	31.5 in	23.5 in	
<b>CM:</b>	80 cm	60 cm	
<b>Individual Product Weight</b>	1.6kg		
<b>Export Packaging:</b>	Corrugated Carton		
<b>Quantity Per Carton</b>	4		
<b>Measurements</b>	<b>L</b>	<b>W</b>	<b>H</b>
<b>Inches:</b>	101.2	17.7	15
<b>CM:</b>	57 cm	45 cm	38 cm
<b>Cube:</b>	0.044cbm		
<b>GW:</b>	8.5Kg		
<b>NW:</b>	7.5kg		
<b>Features:</b>			
<b>Neoprene:</b>	Material: CR Neoprene Density: 0.2 +/- 0.02 Thickness: 4mm		
<b>Jersey:</b>	Material: 100% Nylon		
<b>(Green Outside):</b>			
<b>(Black Inside):</b>	Material: 100% nylon		
<b>Cutting Direction:</b>	Warp		
<b>Printing:</b>	Plastic Ink		
<b>Stitching:</b>	100% Nylon		
<b>Product Testing:</b>	Tensile Testing: (04.73.6002) Laminate Peel Testing: (04.73.6005) Wash Testing and Stretch and Recovery Testing: Concluding the LifeWrap NASG has a lifespan of at least 144 uses		



points/size	small (cm)	large (cm)
panel 6 - tummy total width	110.5	120
panel 4 - hip total width	97.8	96.5
total length - orange line	101	123



**Estimated Shipping & Landed Costs (for ref only)**

Destination:	By Sea via 20' Container			
# of Units	1000 / 20' FCL	500	250	100
Est Freight Cost				
Cost Per Unit				
Estimated Transit / Shipping Time				
Total Cost per Qty Incl Shipping				

Destination:	By Sea via 20' Container			
# of Units	1000 / 20' FCL	500	250	100
Est Freight Cost				
Cost Per Unit				
Estimated Transit / Shipping Time				
Total Cost per Qty Incl Shipping				

Destination:	By Sea via 20' Container			
# of Units	1000 / 20' FCL	500	250	100
Est Freight Cost				
Cost Per Unit				
Estimated Transit / Shipping Time				
Total Cost per Qty Incl Shipping				

