



SALES@VIAGLOBALHEALTH.COM

WWW.VIAGLOBALHEALTH.COM



Innovative

The LifeWrap is a low-cost, low-technology, first-aid device to treat hypovolemic shock



Integrated

Can be used with other interventions, kept on during transportation procedures and surgery




Practical

The suit is affordable, safe, and accessible-ideal for low resource settings


Product Fact Sheet for the LifeWrap NASG

Date:

Description:	LifeWrap NASG Non-Pneumatic Anti-shock Garment		
Size	LARGE / GREEN		
Product Code:	LWNASGLGE15		
Manufacturing MOQ	1,000 Units		
FOB USD\$ FOB HK			
Below Production MOQ	50 Units		
FOB USD\$ HK	POA		
Bar Code #:	0680569839954		
Individual Packaging:	Polybag		
Ind Packed Measurements	L	W	
Inches:	31.5 in	23.5 in	
CM:	80 cm	60 cm	
Individual Product Weight	1.6kg		
Export Packaging:	Corrugated Carton		
Quantity Per Carton	4		
Measurements	L	W	H
Inches:	101.2	17.7	15
CM:	57 cm	45 cm	38 cm
Cube:	0.044cbm		
GW:	8.5Kg		
NW:	7.5kg		
Features:			
Neoprene:	Material: CR Neoprene Density: 0.2 +/- 0.02 Thickness: 4mm		
Jersey: (Green Outside):	Material: 100% Nylon		
(Black Inside):	Material: 100% nylon		
Cutting Direction:	Warp		
Printing:	Plastic Ink		
Stitching:	100% Nylon Tensile Testing: (04.73.6002) Laminate Peel Testing: (04.73.6005)		
Product Testing:	Wash Testing and Stretch and Recovery Testing: Concluding the LifeWrap NASG has a lifespan of at least 144 uses		



points/size	small (cm)	large (cm)
panel 6 - tummy total width	110.5	120
panel 4 - hip total width	97.8	96.5
total length - orange line	101	123



Estimated Shipping & Landed Costs (for ref only)

Destination:	By Sea via 20' Container			
# of Units	1000 / 20' FCL	500	250	100
Est Freight Cost				
Cost Per Unit				
Estimated Transit / Shipping Time				
Total Cost per Qty Incl Shipping				

Destination:	By Sea via 20' Container			
# of Units	1000 / 20' FCL	500	250	100
Est Freight Cost				
Cost Per Unit				
Estimated Transit / Shipping Time				
Total Cost per Qty Incl Shipping				

Destination:	By Sea via 20' Container			
# of Units	1000 / 20' FCL	500	250	100
Est Freight Cost				
Cost Per Unit				
Estimated Transit / Shipping Time				
Total Cost per Qty Incl Shipping				

