



Pumani

bubble CPAP



Product Catalogue

© 2018



Table of Contents

Pumani bubbleCPAP..... Page 3-5

What is included in the package?

What is included in the “Starter Kit”?

Pumani Accessory Items..... Page 6-11

What are the various patient accessories?

How to order patient accessories and spare items?

Product Recommendations..... Page 12-14

Nasal prong and hat size guide

Recommended accessory order quantities

Pumani bubble CPAP



What is included in the package?

The Pumani bubble CPAP package has:

- One Pumani bubble CPAP device
- Pumani 'Starter kit' (see next page)
(Item ID 3-597)



What is included in the 'Starter Kit'?



GLOBAL HEALTH

SOURCE. PURCHASE. SHIP

SALES@VIAGLOBALHEALTH.COM

WWW.VIAGLOBALHEALTH.COM



Patient Tubing

Quantity: 3



Bottle

Quantity: 1



Fish Elbow

Quantity: 6



Hat Clips

Quantity: 6



Hats (Small, Medium, Large)

Quantity: 3

Small - 1 Medium - 1 Large - 1



Variety Nasal Prong Pack

Quantity: 6 prongs, 1 each

Size 0 - 1 Size 3 - 1
Size 1 - 1 Size 4 - 1
Size 2 - 1 Size 5 - 1



Tube Holder

Quantity: 3



Power Cord, UK Style

Quantity: 1



User Manual

Quantity: 1

Pumani Accessory Items



Pumani Patient Accessories

There are 6 unique patient accessories as shown below. Each of these items can be ordered individually.




GLOBAL HEALTH

SOURCE. PURCHASE. SHIP

SALES@VIAGLOBALHEALTH.COM

WWW.VIAGLOBALHEALTH.COM



Nasal Prong

Hat

Tube Holder

Patient Tubing

Hat Clip

Fish Elbow

Patient Accessories - Nasal Prong

Nasal Prongs can be ordered in single size packs or variety size packs.



GLOBAL HEALTH

SOURCE. PURCHASE. SHIP

SALES@VIAGLOBALHEALTH.COM

WWW.VIAGLOBALHEALTH.COM



000 Size 000




Nasal Prong
(Size 000)

Pack of 6

Item ID: 3-548

0 Size 0



Nasal Prong
(Size 0)

Pack of 6

Item ID: 3-520

1 Size 1



Nasal Prong
(Size 1)

Pack of 6

Item ID: 3-521

2 Size 2



Nasal Prong
(Size 2)

Pack of 6

Item ID: 3-522

3 Size 3



Nasal Prong
(Size 3)

Pack of 6

Item ID: 3-523

4 Size 4



Nasal Prong
(Size 4)

Pack of 6

Item ID: 3-524

5 Size 5



Nasal Prong
(Size 5)

Pack of 6

Item ID: 3-526

0-5 Size 0-5
(variety)



Nasal Prong
Variety (Size 0-5)

Pack of 6 (1 each)

Item ID: 3-527

Patient Accessories - Hat



GLOBAL HEALTH

SOURCE. PURCHASE. SHIP

SALES@VIAGLOBALHEALTH.COM

WWW.VIAGLOBALHEALTH.COM

Pumani Accessory Items – Ordering Instructions



Pumani Hat, Size - Small

Pack of 10

Item ID: 3-533



Pumani Hat, Size - Medium

Pack of 10

Item ID: 3-534



Pumani Hat, Size - Large

Pack of 10

Item ID: 3-535

Patient Accessories Continued.,



Hat Clip Pack

Pack of 10

Item ID: 3-536



Tube Holder Pack

Pack of 10

Item ID: 3-552



Fish Elbow Pack

Pack of 10

Item ID: 3-530



Patient Tubing

1 set

Item ID: 3-551



GLOBAL HEALTH

SOURCE. PURCHASE. SHIP

SALES@VIAGLOBALHEALTH.COM

WWW.VIAGLOBALHEALTH.COM

Pumani Spare Items



Bottle

1 piece

Item ID: 3-549



Filter

1 piece

Item ID: 3-541



Flow Meter- Oxygen

1 piece

Item ID: 3-543



Flow Meter-Total

1 piece

Item ID: 3-544



GLOBAL HEALTH

SOURCE. PURCHASE. SHIP

SALES@VIAGLOBALHEALTH.COM

WWW.VIAGLOBALHEALTH.COM

Product Recommendations



Nasal Prong Size and Hat Size Recommendations



GLOBAL HEALTH

SOURCE. PURCHASE. SHIP

SALES@VIAGLOBALHEALTH.COM

WWW.VIAGLOBALHEALTH.COM

Pumani Nasal Prong Size Selection Chart

Patient Weight Range	Nasal Prong Size
Less than 1,000 grams	000 or 0
1,000 grams to 1,250 grams	1
1,250 grams to 2,000 grams	2
2,000 grams to 3,000 grams	3
3,000 grams to 4,000 grams	4
Over 4,000 grams	5

Pumani Hat Size Selection Chart

Patient Weight Range	Hat Size
Less than 1,500 grams	Small
1,500 grams to 3,000 grams	Medium
Over 3,000 grams	Large

Suggested Accessory Order Quantities

Accessory Items	Item #	10 Patients	25 Patients	Requested Quantity
Nasal Prong Pack, Size 000 (Pack of 6)*	3-548	1	1	
Nasal Prong Pack, Size 0 (Pack of 6)*	3-520	1	1	
Nasal Prong Pack, Size 1 (Pack of 6)*	3-521	1	1	
Nasal Prong Pack, Size 2 (Pack of 6)*	3-522	1	1	
Nasal Prong Pack, Size 3 (Pack of 6)*	3-523	1	1	
Nasal Prong Pack, Size 4 (Pack of 6)*	3-524		1	
Nasal Prong Pack, Size 5 (Pack of 6)*	3-526			
Nasal Prong Pack, Variety (Pack of 6)*	3-527	1	1	
Small Hat (Pack of 10)**	3-533	1	2	
Medium Hat (Pack of 10)**	3-534	1	2	
Large Hat (Pack of 10)**	3-535	1	1	
Hat Clip (Pack of 10)	3-536	1	1	
Patient Tubing*	3-551	10	25	
Tube Holder (Pack of 10)	3-552	1	2	
Fish Elbow (Pack of 10)*	3-530	2	5	
Bottle^	3-549	1	1	
Filter^	3-541	1	1	

Suggested quantities are based on the estimated number of patients to be treated per year and may need to be adjusted for individual hospitals.

*Items labeled as single-use only
 **Hats may need to be changed every couple days to ensure they remain secure on the head
 ^Replacement parts recommended yearly per Pumani CPAP device



GLOBAL HEALTH

SOURCE. PURCHASE. SHIP

SALES@VIAGLOBALHEALTH.COM

WWW.VIAGLOBALHEALTH.COM



Pumani
bubble CPAP

Product Catalogue

© 2018



GLOBAL HEALTH

SOURCE. PURCHASE. SHIP

SALES@VIAGLOBALHEALTH.COM

WWW.VIAGLOBALHEALTH.COM